



Code	Immunoturbidimetry - IMT	Packaging	Specifications
IMT279-080	α1–ACID GLYCOPROTEIN	R1 3x25 ml R2 1x5 ml	4 mg/dl - 300 mg/dl
IMT276-080	Antitrombin III	R1 3x25 ml R2 1x5 ml	8 mg/dl - 65 mg/dl
IMT481-050	ASO ANTI-STREPTOLYSIN O (with calibrator)	R1 2x20 ml R2 1x10 ml CAL 1x2 ml	3 IU/ml - 800 IU/ml
IMT581-100	ASO ANTI-STREPTOLYSIN O (with calibrator)	R1 4x20 ml R2 2x10 ml CAL 1x2 ml	3 IU/ml - 800 IU/ml
IMT277-065	C 3	R1 3x20 ml R2 1x5 ml	20 mg/dl - 400 mg/dl
IMT278-065	C 4	R1 3x20 ml R2 1x5 ml	2 mg/dl - 80 mg/dl
IMT483-050	CRP C-REACTIVE PROTEIN (with calibrator)	R1 2x20 ml R2 1x10 ml CAL 1x1ml	1.0 mg/l - 150 mg/l
IMT583-100	CRP C-REACTIVE PROTEIN (with calibrator)	R1 4x20 ml R2 2x10 ml CAL 1x1ml	1.0 mg/l - 150 mg/l
IMT458-036	Direct HbA1c (130 test)	R1 1x30 ml R2 1x6 ml	16 g/dl
IMT584-045	FERRITIN (with calibrator)	R1 2x15 ml R2 1x15 ml CAL 1x3ml	4 µg/l - 500 µg/l
IMT289-080	FIBRINOGEN	R1 3x25 ml R2 1x5 ml	40 mg/dl - 500 mg/dl
IMT282-045	IgA	R1 2x20 ml R2 1x5 ml	1 mg/dl - 600 mg/dl
IMT281-055	IgG	R1 2x25 ml R2 1x5 ml	60 mg/dl - 2700 mg/dl
IMT283-045	IgM	R1 2x20 ml R2 1x5 ml	6 mg/dl - 500 mg/dl
IMT286-080	KAPPA LIGHT CHAIN	R1 3x25 ml R2 1x5 ml	40 mg/dl - 800 mg/dl
IMT287-080	LAMBDA LIGHT CHAIN	R1 3x25 ml R2 1x5 ml	20 mg/dl - 450 mg/dl
IMT275-080	LIPASE (UV turbidimetric)	R1a 1x70 ml R2b 20x2,5 ml CAL 4x 3 ml	700 U/l
IMT284-045	MICROALBUMIN (with calibrator)	R1 2x20 ml R2 1x5 ml CAL 1x1 ml	0.7 mg/l - 400 mg/l
IMT482-050	RF RHEUMATOID FACTORS (with calibrator)	R1 2x20 ml R2 1x10 ml CAL 1x2 ml	2 IU/ml - 160 IU/ml
IMT582-100	RF RHEUMATOID FACTORS (with calibrator)	R1 4x20 ml R2 2x10 ml CAL 1x2 ml	2 IU/ml - 160 IU/ml
IMT285-080	TRANSFERRIN	R1 3x25 ml R2 1x5 ml	40 mg/dl - 500 mg/dl

Code	Standards & Calibrators & Controls	Packaging	Specifications
IMTSTD780	FIBRINOGEN STANDARD	1x0,5 ml	The value for Fibrinogen is indicated on label
IMTSTD784-001	PEDIATRIC STANDARD	1x1 ml	Calibrations values: AMG-TRF-IgA-IgG-IgM-KAP-LAM
IMTSTD794-001	PROTEIN STANDARD HIGH	1x1 ml	Calibration values: IgG-IgA-IgM-C1-CER-AT3-AGP-ALB-C3-C4-AAT-HAP-TRF-KAP-LAM-PAL-AMG
IMTSTD796-005	PROTEIN STANDARD SET (5 levels)	5x1 ml	Calibration values: IgG-IgA-IgM-C1-CER-AT3-AGP-ALB-C3-C4-AAT-HAP-TRF-KAP-LAM-PAL-AMG
IMTCAL137-002	HbA1c CALIBRATOR SET	1x2 ml + 4x0,5 ml	Liquid calibrator
IMTCTR128-001	HbA1c CONTROL SET	2x0,5 ml	Set (2 level)
IMTCTR791-001	FIBRINOGEN CONTROL	1x1 ml	250 mg/dl
IMTCTR805-001	IMMUNOLOGY CONTROL LOW	1x1 ml	Control values: ALB-AGP-AAT-AMG-ASL-AT3-BMG-CER-C3-C4-C1-CRP-FRT-HAP-IgA-IgM-IgG-KAP-PAL-RF-TRF
IMTCTR806-001	IMMUNOLOGY CONTROL HIGH	1x1 ml	Control values: ALB-AGP-AAT-AMG-ASL-AT3-BMG-CER-C3-C4-C1-CRP-FRT-HAP-IgA-IgM-IgG-KAP-PAL-RF-TRF
IMTCRT790-001	MICROALBUMIN CONTROL	1x1 ml	200 mg/dl
IMTCTR798-001	RHEUMATHOID CONTROL SERUM LEV 1	1x1 ml	Control values: ASO-CRP-RF
IMTCTR799-001	RHEUMATHOID CONTROL SERUM LEV 2	1x1 ml	Control values: ASO-CRP-RF